Product Data Sheet

AccuPack System

Product Images

Description

FEATURES

• The Accupak system provides a wireless system allowing the user an average of 100 cycles of chair movements. This means the chair can be positioned and used in any area of a room without relying on a nearby power source. When the battery requires charging a beeping sound can be heard to warn the user when any button on the handset is pressed. The transformer can then be plugged into the charging outlet and re-charged for a period of 24 hours for maximum charge. The system is fitted to the base of the chair (out of sight), and the side of the chair is fitted with an outlet for charging.

Bristol Maid: Manufacturer and distributor of medical furniture since 1953. An approved NHS and public sector supplier, our equipment can be found in hospitals throughout the UK and internationally.

Code	Description	Weight
AccuPack System	Allows mains free chair movements	-