

High Back Patient Armchairs - Alderbury

Product Images





Description

FEATURES

- BS EN 16139: 2013 Level 2. Furniture. Strength, durability & safety. Requirements for non-domestic seating
- Ergonomic design to ensure correct sit & posture
- Handgrips assist positioning from standing & seated positions
- Solid beech frame (PEFC accredited)
- Mortice & tenon joint construction, glued & pinned
- Frame finished with anti-bacterial varnish
- Plastic glides fitted to chair legs
- Contoured seat & back
 - o Fixed Seat is fixed to the frame
 - Removable Seat can be removed to assist cleaning (smooth vinyl upholstered underside) & is reversible to allow the seat section to be repositioned to provide additional clearance at the front of the chair
- Hygiene gap prevents the trapping of dirt & bacteria between the seat & back, provides access for cleaning
- Fire retardant foam (BS 5852: Part 2: 1982)

- Upholstered with fire retardant & antimicrobial vinyl
- Maximum capacity 160kg

OVERALL DIMENSIONS (W X D X H)

- Overall Dimensions 660 x 650 x 1195mm
- Seat 450 x 435 / 475* at Height 500mm
- Width Between Arms 530mm
- Arm Height 655mm

STANDARD COLOURS

• Vinyl - Delft

Code	Description	Weight
5ND/HAF/ALDE/VVDEL	Alderbury, High Back Armchair, Fixed Seat, Vinyl Delft	14.4kg
5ND/HAF/ALDE/VVTUR	Alderbury, High Back Armchair, Fixed Seat, Vinyl Türkis	14.4kg
5ND/HAF/ALDE/VVROT	Alderbury, High Back Armchair, Fixed Seat, Vinyl Rot	14.4kg
5ND/HAF/ALDE/VVLIC	Alderbury, High Back Armchair, Fixed Seat, Vinyl Lichtgrau	14.4kg
5ND/HAR/ALDE/VVDEL	Alderbury, High Back Armchair, Removable Seat, Vinyl Delft	14.4kg
5ND/HAR/ALDE/VVTUR	Alderbury, High Back Armchair, Removable Seat, Vinyl Türkis	14.4kg
5ND/HAR/ALDE/VVROT	Alderbury, High Back Armchair, Removable Seat, Vinyl Rot	14.4kg
5ND/HAR/ALDE/VVLIC	Alderbury, High Back Armchair, Removable Seat, Vinyl Lichtgrau	14.4kg

^{*}Seat depth can be altered by repositioning the removable seat pad